## Plástico Essix ACE®

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
## Essix ACE® Plastic

(Instructions for use on Essix® Vacuum Machines only)

### Indications for Use:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

### Intended Use:

Essix ACE® Plastic is thermoformed to fabricate intra-oral appliances such as retainers, aligners for minor tooth movement, and temporary bridges. Essix ACE® Plastic is suitable for use with non-bruxers.

**Contraindications:** None known.

### Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

### Precautions:

Essix ACE® Plastic is not recommended for bruxers.

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only.

**Adverse Reactions:** None known.

### Instructions for Use:

1. Essix ACE® Plastic has a protective film on both sides of the sheet. REMOVE FILM FROM BOTH SIDES BEFORE THERMOFORMING.
2. All models should be trimmed to ¾" high. Full arch models should have the palate trimmed away.
3. **IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
4. Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
5. Place the plastic into the frame of the vacuum machine and the model on the base plate.
6. After pre-heating, swing the heating element over the plastic.
7. Observe the plastic as it heats.
8. Essix ACE® Plastic should be heated until it sags approximately ¾" – ½".
9. Move the heating element away from the plastic.
10. Turn on the vacuum and lower the plastic over the model.
11. Push the plastic into each undercut with the Essix Accentuator® Tool.
12. Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it “shrinks” the material around the cast, realizing a tighter fit.
13. Trim and finish the appliance as desired.

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic “folds” when thermoformed, reduce heating time until material does not fold.

Thickness	Vacuum Heat Time*	Biostar™ Code*	Drufomat™
0.75mm (.030")	25 seconds	113	55 seconds
0.88mm (.035")	30 seconds	123	60 seconds
1.0mm (.040")	35 seconds	133	65 seconds

\*after preheating element

Manufactured by:  
**DENTSPLY International Inc.**  
**DENTSPLY Raintree Essix Glenroe**  
1912 44th Avenue East  
Bradenton, FL 34203  
Tel: 800.883.8733  
www.essix.com

Emergo Europe  
Molenstraat 15,  
2513 BH, The Hague  
The Netherlands  
Tel: (+31) 70 345 8570  
Fax: (+31) 70 346 7299

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