Plástico Essix A+®

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Essix A+® Plastic

(Instructions for use on Essix® Vacuum Machines only)

Indications for Use:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

Intended Use:

Essix A+® Plastic is thermoformed to fabricate intra-oral appliances such as bite planes, splints, and snoring appliances. Essix A+® Plastic is suitable for use with non-bruxers.

Contraindications: None known.

Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

Precautions:

Essix A+® Plastic is not recommended for bruxers.

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only.

Adverse Reactions: None known.

Instructions for Use:

- Essix A+® Plastic has a protective film on both sides. REMOVE FILM FROM BOTH SIDES BEFORE THERMOFORMING.
- All models should be trimmed to ¾" high. Full arch models should have the palate trimmed away.
- IMPORTANT: Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
- Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
- Place the plastic into the frame of the vacuum machine and the model on the base plate.
- After pre-heating, swing the heating element over the plastic.
- Observe the plastic as it heats.
- Essix A+® Plastic should be heated until it sags approximately ½".
- Move the heating element away from the plastic.
- Push the plastic into each undercut with the Essix Accentuator® Tool.
- Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it "shrinks" the material around the cast, realizing a tighter fit.
- Turn on the vacuum and lower the plastic over the model.
- Let vacuum continue to run until plastic is cool.
- Trim and finish the appliance as desired.

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic "folds" when thermoformed, reduce heating time until material does not fold.

| Thickness | Vacuum Heat Time* | Biostar™ Code* | Drufomat™ |
|----------------|-------------------|----------------|-------------|
| 0.5mm (.020") | 30 seconds | 122 | 60 seconds |
| 0.75mm (.030") | 35 seconds | 132 | 60 seconds |
| 1.0mm (.040") | 40 seconds | 142 | 65 seconds |
| 1.5mm (.060") | 50 seconds | 164 | 70 seconds |
| 2.0mm (.080") | 70 seconds | 205 | 90 seconds |
| 3.0mm (.120") | 80 seconds | 226 | 130 seconds |

*after preheating element

Manufactured by:
DENTSPLY International Inc.
DENTSPLY Raintree Essix Glenroe
 1912 44th Avenue East
 Bradenton, FL 34203
 Tel: 800.883.8733
 www.essix.com

EC REP Emergo Europe
 Molenstraat 15,
 2513 BH, The Hague
 The Netherlands
 Tel: (+31) 70 345 8570
 Fax: (+31) 70 346 7299

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