Essix Embrace® Plastic

(Instructions for use on Essix® Vacuum Machines Only)

Indications for Use:

Essix Embrace® Plastic is thermoformed to fabricate intra-oral appliances such as bite planes, splints, snoring appliances and nightguards. Essix Embrace® Plastic is not suitable for use with bruxers.

Contraindications: None known.

Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

Precautions:

Not recommended for bruxers.

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in cool, dry place.

Single Use Only.

Adverse Reactions: None known.

Directions for Use:

1. All models should be trimmed to ¾" high. Full arch models should have the palate trimmed away.
2. IMPORTANT: Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
3. Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
4. Place the Essix Embrace® plastic into the frame of the vacuum machine and the model on the base plate. NOTE: Place curved side up into the frame.
5. After pre-heating, swing the heating element over the plastic.
6. Observe the plastic as it heats.
7. Essix Embrace® will completely flatten – this is the optimal thermoforming time.
8. Move the heating element away from the plastic.
9. Turn on the vacuum and lower the plastic over the model.
10. Push the plastic into each undercut with the Essix Accentuator® Tool.
11. Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it “shrinks” the material around the cast, realizing a tighter fit.
12. Let vacuum continue to run until plastic is cool.
13. Trim and finish appliance as desired.

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic “folds” when thermoformed, reduce heating time until material does not fold.

Thickness	Vacuum Heat Time*	Biostar™ Code*	Drufomat™
0.75mm (.030")	35 seconds	133	70 seconds
1.0mm (.040")	40 seconds	143	75 seconds

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