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**DENTSPLY International Inc.**  
**DENTSPLY Raintree Essix Glenroe**

1912 44th Avenue East  
Bradenton, FL 34203  
Tel: 800.883.8733  
www.essix.com

Biostar y Drufomat no son marcas comerciales de DENTSPLY.

ESP #MKNL RTE-193-11 Rev. 09/11

## Essix® Nightguard Laminate Plastic

(Instructions for use on Essix® Vacuum Machines only)

### Indications for Use:

Mouthguard and aligner materials are indicated for the fabrication of orthodontic and dental appliances such as aligners, bite planes, mouthguards, nightguards, snoring appliances, splints, retainers, repositioners, and temporary bridges.

### Intended Use:

Essix® Nightguard Laminate Plastic is thermoformed to fabricate intra-oral appliances such as bruxism appliances, nightguards, orthodontic repositioners, and periodontal splints. Essix® Nightguard Laminate Plastic is suitable for use with bruxers/grinders.

**Contraindications:** None known.

### Warnings:

Excessive heating time during thermoforming will cause plastic to over sag and potential for injury to operator.

Improper thermoforming of the plastic appliance may cause the appliance to crack/break, resulting in sharp edges, loose pieces and possible aspiration of pieces.

### Precautions:

Patients with a history of allergic reactions to plastics should not use this product.

Store plastic sheets in a cool, dry place.

Single Use Only.

**Adverse Reactions:** None known.

### Directions for Use:

1. All models should be trimmed to ¾" high. Full arch models should have the palate trimmed away.
2. **IMPORTANT:** Pre-heat the machine before using! Pre-heat the vacuum machine for at least 3 minutes or until the heating element is a solid orange color. (Heating times may vary with current fluctuations or length of time the heating element has been left on.)
3. Spray the model with Trim-Rite® Release Agent to prevent the material from sticking to the model.
4. Place the plastic into the frame on the vacuum machine and the model on the base plate.
5. Swing the heating element over the plastic. Observe the material as it heats.
6. Material will sag approximately ½" – this is the optimal thermoforming time.
7. Move the heating element away from the plastic.
8. Turn on the vacuum and lower the plastic over the model.
9. Immediately spray the plastic encased model with Essix Freeze Spray® Coolant. Not only does it cool the plastic, prohibiting further thinning, but it "shrinks" the material around the cast, realizing a tighter fit.
10. Trim and finish appliance as desired.

Heating time may vary depending upon machine. If plastic does not adapt to model, add heating time until adaptation occurs. If plastic "folds" when thermoformed, reduce heating time until material does not fold.

Thickness	Vacuum Heat Time*	Biostar™ Code*	Drufomat™
3.0mm (.120")	50 seconds	225	120 seconds

\*after preheating element

Manufactured by:

**DENTSPLY International Inc.**  
**DENTSPLY Raintree Essix Glenroe**

1912 44th Avenue East  
Bradenton, FL 34203  
Tel: 800.883.8733  
www.essix.com

Biostar and Drufomat are not trademarks of DENTSPLY.

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